

MiniLap System with MiniGrip Handle

Smaller by Design

Teleflex is committed to helping make minimally invasive surgery less invasive through smaller incisions for less trauma and smaller scars.

Introducing the MiniLap Percutaneous Surgical System with MiniGrip Handle from Teleflex:

Designed by a surgeon, the slim 2.4 mm shaft diameter instruments can be percutaneously inserted using an integrated needle tip. The grasper jaws can then be deployed to grasp and manipulate tissue. The patented deployment design allows grasper jaws to open up to 12.5 mm to grasp difficult structures, like distended gall bladders, and the proprietary steel shaft provides the strength and rigidity needed for secure tissue retraction.

MiniLap can be used for a variety of procedures such as lap cholecystectomy and lap appendectomy, and benign gynecology procedures such as diagnostic laparoscopy, LSO, and endometriosis.


Fewer Trocars

- The integrated needle tip allows direct insertion through the skin, meaning fewer trocars are used.
- In many cases, 2 or more trocars can be eliminated.

Less Trauma

- The sturdy 2.4 mm shaft allows for insertion directly through the skin, leaving a nearly unnoticeable scar.
- With MiniLap, you can improve patient satisfaction while achieving the outcomes you want for your patients.

Diameter Comparison



Benefits



The Clinician
Less trauma. Easy to Use.
Sturdy.




Your Institution
Fewer trocars. Patient satisfaction.



The Patient
Nearly unnoticeable scars and potentially less pain.



MiniLap System with MiniGrip Handle



2.4 mm Shaft allows for insertion directly through the skin, leaving a minimal scar

Jaws open up to 12.5 mm aperture

MiniLap Percutaneous Graspers



2.3 mm Shaft allows for insertion directly through the skin, leaving a minimal scar

Jaws open up to 12.5 mm aperture

MiniPolar Electrosurgery Probes



2.4 mm Shaft allows for insertion directly through the skin, leaving a minimal scar

MiniLap System with MiniGrip Handle

ITEM NUMBER	PRODUCT DESCRIPTION	UNIT PER BOX
PGAC300	MiniGrip Alligator Grasper	6
PGCC300	MiniGrip Clutch Grasper	6
PGBK300	MiniGrip Babcock Grasper	6
PGBC300	MiniGrip Tong Grasper	6

MiniLap Percutaneous Surgical System

ITEM NUMBER	PRODUCT DESCRIPTION	UNIT PER BOX
GBC250	Alligator Grasper	6
CLC250	Clutch Grasper	6
BCK250	Babcock Grasper	6
BLC250	Bowel Grasper	6

MiniPolar Percutaneous Electrosurgical Probes

ITEM NUMBER	PRODUCT DESCRIPTION	UNIT PER BOX
ECMC300	Curved Spatula Probe	6
ECMS300	Straight Spatula Probe	6
ECMH300	Hook Probe	6
ECMP300	Conical Probe	6

To order:

Please contact us at orders.intl@teleflex.com

Teleflex, the Teleflex logo, MiniGrip, MiniLap and MiniPolar are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this material is not a substitute for the product Instructions for Use. Not all products may be available in all countries. Please contact your local representative. Revised: 08/2017. © 2017 Teleflex Incorporated. All rights reserved. 94 08 17 - 00 00 01 - REV A - 08 17 02.5